
UM-SJTU JOINT INSTITUTE
PHYSICS LABORATORY
(VP241)

LABORATORY REPORT

EXERCISE 4
POLARIZATION OF LIGHT

Name: **Kang Jiaming**
Name: Hong Yichen

ID: **518021911220**
ID: 518370910011

Group: 17
Group: 17

Date: 8 Nov. 2019

Contents

| | | |
|----------|--|-----------|
| 1 | Introduction | 1 |
| 1.1 | Polarization of Light | 1 |
| 1.2 | Polarizer | 1 |
| 1.3 | Malus' Law | 1 |
| 1.4 | Generation of Elliptically and Circularly Polarized Light. Half-wave and Quarter-wave Plates | 2 |
| 2 | Experimental Setup | 4 |
| 2.1 | Apparatus | 4 |
| 3 | Measurement Procedure | 4 |
| 3.1 | Apparatus Adjustment | 4 |
| 3.2 | Demonstration of Malus' Law | 4 |
| 3.3 | Linearly Polarized Light and the Half-wave Plate | 4 |
| 3.4 | Circularly and Elliptically Polarized Light and the 1/4-wave Plate | 5 |
| 4 | Results | 5 |
| 4.1 | Demonstration of Malus' Law | 5 |
| 4.2 | Linearly Polarized Light and the Half-wave Plate | 7 |
| 4.3 | Circularly and Elliptically Polarized Light and the 1/4-wave Plate | 8 |
| 4.3.1 | Rotation Angle: 0° | 8 |
| 4.3.2 | Rotation Angle: 20° | 11 |
| 4.3.3 | Rotation Angle: 45° | 14 |
| 4.3.4 | Rotation Angle: 70° | 18 |
| 5 | Discussion and Conclusions | 18 |
| 5.1 | Demonstration of Malus' Law | 18 |
| 5.2 | Linearly Polarized Light and the Half-wave Plate | 19 |
| 5.3 | Circularly and Elliptically Polarized Light and the 1/4-wave Plate | 19 |
| 5.3.1 | Rotation Angle: 0° | 19 |
| 5.3.2 | Rotation Angle: 20° | 19 |
| 5.3.3 | Rotation Angle: 45° | 19 |
| 5.3.4 | Rotation Angle: 70° | 19 |
| 5.4 | Causes for errors and uncertainties | 19 |
| 5.5 | Suggestions | 20 |
| 5.6 | Conclusions | 20 |
| A | Measurement Uncertainty Analysis | 21 |
| A.1 | Uncertainty for Data in Demonstration of Malus' Law | 21 |
| A.2 | Uncertainty for Linearly Polarized Light and the Half-wave Plate | 21 |
| A.3 | Uncertainty for Circularly and Elliptically Polarized Light and the 1/4-wave Plate | 22 |
| B | Data Sheet | 23 |

1 Introduction

The objective of this exercise is to study the polarization of light. Specifically, Malus' law will be verified and the way how half- and quarter-wave plates work in optical systems will be explored.

1.1 Polarization of Light

Light in nature is electromagnetic wave. The electric field vector \mathbf{E} in the electromagnetic wave is referred to as the light vector. In the plane perpendicular to the propagation direction of a light wave, the light vector may have different directions along which its magnitude oscillates. If the light vector oscillates in all possible transverse direction, the light is called *natural light*. Otherwise, the light is called *polarized*. The light, for which the light vector maintains a certain oscillation direction, is called *linearly polarized* and the axis defining the direction is called the polarization axis. The light with the light vector direction rotating about the propagation direction, so that its endpoint traces a circle, is called *circularly polarized* light. If the vector traces an ellipse, the light is said to be *elliptically polarized*.

1.2 Polarizer

Polarizer is a device that can produce polarized light. It polarizes the light by only allowing the light polarized in a certain direction to pass through. With such a property, it can also be used to detect and analyze the polarization state of light (it is then called an analyzer).

1.3 Malus' Law

A fundamental law that the polarization of light obeys is as follows. Suppose that we have two polarizers arranged so that their planes are parallel — the left one plays the role of a polarizer, the other one is an analyzer (see Figure 1). Let the angle between their transmission directions (polarization axes) be θ . The light is incident normally on the polarizer and then continues to the analyzer. The intensity of the linearly polarized light leaving the analyzer is

$$I_{\text{light}} = I_{\text{light},0} \cos^2 \theta, \quad (1)$$

where $I_{\text{light},0}$ is the intensity of the linearly polarized light incident on the analyzer. This is the so-called Malus' law.




Figure 1: Change in the brightness of the light depends on the mutual orientation of the polarizer and the analyzer.

1.4 Generation of Elliptically and Circularly Polarized Light. Half-wave and Quarter-wave Plates

***This section is selected from [1], but the crucial points regarding this lab are highlighted here.**

Suppose that a linearly polarized light is incident normally on a crystal plate whose surface is parallel to its optical axis, and the angle between the polarizing axis and the optical axis of the plate is α . Then the linearly polarized light is resolved into two waves: an e-wave with the oscillation direction parallel to the optical axis of the plate (extraordinary axis) and an o-wave whose oscillation direction is perpendicular to the optical axis (ordinary axis). They propagate in the same direction, but with different speeds. The resulting optical path difference over the thickness d of the plate is

$$\Delta = (n_e - n_o)d,$$

and, consequently, the phase difference

$$\delta = \frac{2\pi}{\lambda}(n_e - n_o)d,$$

where λ is the wavelength, n_e is the refractive index for the extraordinary axis, and n_o is the refractive index for the ordinary axis. In a so-called positive crystal $\delta > 0$, whereas in a negative one $\delta < 0$.

As shown in Figure 4, when the light propagates through the crystal plate, the two components of the light vector are

$$\begin{aligned} E_x &= A_o \cos \omega t \\ E_y &= A_e \cos(\omega t + \delta), \end{aligned}$$

where $A_e = A \cos \alpha$, $A_o = A \sin \alpha$. Eliminating time from the above equations one obtains

$$\frac{E_x^2}{A_o^2} + \frac{E_y^2}{A_e^2} - 2 \frac{E_x E_y}{A_o A_e} \cos \delta = \sin^2 \delta, \quad (2)$$

which is the equation of an ellipse.




Figure 2: Linearly polarized light passing through a waveplate.

When the thickness of the plate changes, the optical path difference changes as well. Some cases of particular interest, are discussed below:

- If $\Delta = k\lambda$, where $k = 0, 1, 2, \dots$, the phase difference $\delta = 0$, and Eq. (2) reduces to

$$E_y = \frac{A_e}{A_o} E_x,$$

which is a linear equation. Hence the transmitted light is linearly polarized with the oscillation direction remaining unchanged. A waveplate that satisfies this condition is called a *full-wave plate*. The light goes through a full-wave plate without changing its polarization state.

- If $\Delta = (2k + 1)\lambda/2$, where $k = 0, 1, 2, \dots$, the phase difference $\delta = \pi$, and Eq. (2) simplifies to

$$E_y = -\frac{A_e}{A_o} E_x.$$

The transmitted light is also linearly polarized with the polarization axis rotated by the angle of 2α . A waveplate that satisfies the condition is called *1/2-wave plate* or *half-wave plate*. When a polarized light passes through a half-wave plate, its polarization axis gets rotated by an angle 2α . If $\alpha = \pi/4$, then the polarization axis of the transmitted light is perpendicular to that of the incident light.

- Finally, if $\Delta = (2k + 1)\lambda/4$, where $k = 0, 1, 2, \dots$, the phase difference $\delta = \pm\pi/2$, and Eq. (2) transforms into

$$\frac{E_x^2}{A_o^2} \pm \frac{E_y^2}{A_e^2} = 1.$$

The transmitted light is elliptically polarized with a waveplate that satisfies the above condition is called a *1/4-wave plate* or a *quarter-waveplate* and is an important optical element in many polarization experiments.

If $A_e = A_o = A$, then $E_x^2 + E_y^2 = A^2$, and the transmitted light is circularly polarized. Since the amplitudes of the *o*-wave and the *e*-wave are both functions of α , the polarization state after passing through a 1/4-wave plate will vary, depending on the angle:

- if $\alpha = 0$, the transmitted light is linearly polarized with the polarization axis parallel to the optical axis of the 1/4-wave plate;

- if $\alpha = \pi/2$, the transmitted light is linearly polarized with the polarization axis perpendicular to the optical axis of the 1/4-wave plate;
- if $\alpha = \pi/4$, the transmitted light is circularly polarized;
- otherwise, the transmitted light is elliptically polarized.

2 Experimental Setup

2.1 Apparatus

The experimental setup consists of a laser, a silicon photo-cell, a UT51 digital universal meter, two polarizers, a 1/2-wave plate, a 1/4-wave plate and an optical bench where the elements are placed.

The precision of the devices is shown in Table 1.

| Instrument | Quantity | Precision |
|----------------------|----------------|-----------------------|
| Scale on the element | Angle θ | 2° |
| Universal meter | Current I | $0.001 [\mu\text{A}]$ |

Table 1: Precision of the measurement instruments.

3 Measurement Procedure

3.1 Apparatus Adjustment

1. Adjust the laser and the photo-cell so that the light can pass through the ϕ 6.0 aperture.
2. Placed the analyzer onto the optical bench, adjust its position so that the current detected is in the range $0.8 \sim 1.8 \mu\text{A}$.

3.2 Demonstration of Malus' Law

1. Place an analyzer between the polarizer and the photo-cell. Record this value as I_0 .
2. Adjust the angle of the analyzer until the electric current reaches 0. Set this position as $\theta = 90^\circ$.
3. Rotate the analyzer from 90° to 0° and record the magnitude of the current I every 5° and record the values.
4. To verify Malus' Law, plot the graph I/I_0 vs. $\cos^2 \theta$ and perform linear fitting.

3.3 Linearly Polarized Light and the Half-wave Plate

1. Place the 1/2-wave plate between the polarizer and the analyzer. Rotate it to make the light extinction appear again and set this position as the initial position.
2. Rotate the 1/2-wave plate for $\alpha = 10^\circ$ from the initial position and rotate the analyzer to make the light extinction appear again, record the angle of rotation $\Delta\theta$.
3. To explore the relationship, plot the graph $\Delta\theta$ vs. θ .

3.4 Circularly and Elliptically Polarized Light and the 1/4-wave Plate

1. Replace the 1/2-wave plate by the 1/4-wave plate and rotate it to make the light extinction appear and set this position as the initial position. At this point $\alpha = 0^\circ$. Then rotate the 1/4-wave plate and observe the change in the light intensity.
2. Rotate the analyzer for 360° and record the light intensity (which is indicated by the current I) for every 10° and record the data in a table.
3. Rotate the 1/4-wave plate for 20° , repeat Step 3.
4. Rotate the 1/4-wave plate for 45° , and repeat Step 3.
5. Rotate the 1/4-wave plate for 70° . Rotate the analyzer and record its position and the magnitude of the current when the light intensity reaches a maximum.
6. To find out the relation, plot the relation between the rotation angle of the analyzer and the light amplitude in polar coordinates. Normalize the amplitude by its maximum value. Mark the position recorded in Step 6 and compare it with the data recorded.
7. To compare with the circular polarization, plot a linear fit to the data when the angle is 45° .

Cautions

1. Do not direct the laser beam into the eye.
2. Do not touch the surface of the polarizers or the wave plates.

4 Results

4.1 Demonstration of Malus' Law

The measurement data are presented in Table 2.

| Maximum Electric Current I_0 | | $1.037 \pm 0.001 [\mu\text{A}]$ | |
|---------------------------------|---|---------------------------------|---|
| $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ | $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ |
| 0 | 1.037 | 50 | 0.470 |
| 5 | 1.033 | 55 | 0.373 |
| 10 | 1.014 | 60 | 0.294 |
| 15 | 0.986 | 65 | 0.206 |
| 20 | 0.947 | 70 | 0.147 |
| 25 | 0.896 | 75 | 0.085 |
| 30 | 0.822 | 80 | 0.042 |
| 35 | 0.730 | 85 | 0.015 |
| 40 | 0.643 | 90 | 0.003 |
| 45 | 0.560 | | |

Table 2: Measurement data for Malus' law demonstration.

To find the relation between $\cos^2 \theta$ and I/I_0 , the two set of values need to be calculated first. Take the first set of data as an example.

$$\cos^2 \theta = \cos^2(0 \times \frac{\pi}{180}) = 1 \pm 0,$$

$$\frac{I}{I_0} = \frac{1.037}{1.037} = 1.0000 \pm 0.0013.$$

Perform similar calculations to each of the other sets of data and the results are shown in Table 3.

| $\theta [^\circ] \pm 2 [^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ | $\cos^2 \theta$ | I/I_0 |
|----------------------------------|---|-------------------|---------------------|
| 0 | 1.037 | 1 ± 0 | 1.0000 ± 0.0013 |
| 5 | 1.033 | 0.992 ± 0.006 | 0.9961 ± 0.0013 |
| 10 | 1.014 | 0.970 ± 0.012 | 0.9778 ± 0.0013 |
| 15 | 0.986 | 0.933 ± 0.017 | 0.9508 ± 0.0013 |
| 20 | 0.947 | 0.88 ± 0.02 | 0.9132 ± 0.0013 |
| 25 | 0.896 | 0.82 ± 0.03 | 0.8640 ± 0.0013 |
| 30 | 0.822 | 0.75 ± 0.03 | 0.7927 ± 0.0012 |
| 35 | 0.730 | 0.67 ± 0.03 | 0.7040 ± 0.0012 |
| 40 | 0.643 | 0.59 ± 0.03 | 0.6201 ± 0.0011 |
| 45 | 0.560 | 0.50 ± 0.03 | 0.5400 ± 0.0011 |
| 50 | 0.470 | 0.41 ± 0.03 | 0.4532 ± 0.0011 |
| 55 | 0.373 | 0.33 ± 0.03 | 0.3597 ± 0.0010 |
| 60 | 0.294 | 0.25 ± 0.03 | 0.2835 ± 0.0010 |
| 65 | 0.206 | 0.18 ± 0.03 | 0.1986 ± 0.0010 |
| 70 | 0.147 | 0.12 ± 0.02 | 0.1418 ± 0.0010 |
| 75 | 0.085 | 0.067 ± 0.017 | 0.0820 ± 0.0010 |
| 80 | 0.042 | 0.030 ± 0.012 | 0.0405 ± 0.0010 |
| 85 | 0.015 | 0.008 ± 0.006 | 0.0145 ± 0.0010 |
| 90 | 0.003 | 0 ± 0 | 0.0029 ± 0.0010 |

Table 3: Result for $\cos^2 \theta$ and I/I_0 .

Then, a linear fit of the plot I/I_0 vs. θ are performed (Figure 3). The slope of the linear fitting is 1.01 ± 0.02 .




Figure 3: Linear fit of I/I_0 vs. $\cos^2 \theta$ relation.

4.2 Linearly Polarized Light and the Half-wave Plate

The measurement data are presented in Table 4. Note that we have substracted each of the data by the initial angle in the report.

| Rotation angle of the 1/2-wave plate [$^\circ$] $\pm 2[^\circ]$ | Rotation angle of the analyzer [$^\circ$] $\pm 2[^\circ]$ |
|---|---|
| initial | 0 |
| 10 | 22 |
| 20 | 40 |
| 30 | 61 |
| 40 | 82 |
| 50 | 102 |
| 60 | 121 |
| 70 | 142 |
| 80 | 162 |
| 90 | 181 |

Table 4: Measurement data for the 1/2-wave plate.

To find the relation between $\Delta\theta$ and θ , the data are plotted in Figure 4 and linear fit is performed. The slope of the linear fit is 2.01 ± 0.02 .




Figure 4: Linear fit of $\Delta\theta$ vs. θ .

Besides, in this part, when the half-wave plate rotates for 360° , 4 times of light extinction are observed. And when the analyzer rotates 360° , 2 times of light extinction during the analyzer rotating 360° . It can be concluded from the above result that **the polarization axis get rotated by twice of the origin angle (2α)** after the light passes through the 1/2-wave plate,

4.3 Circularly and Elliptically Polarized Light and the 1/4-wave Plate

4.3.1 Rotation Angle: 0°

The measurement data for 0° rotation angle of 1/4-wave plate are presented in Table 5. Note that when filling the data sheet, 0° and 90° are mistaken and the mistake is corrected in the report.

| Rotation angle of 1/4-wave plate: 0° | | | |
|--------------------------------------|---|---------------------------------|---|
| Maximum Electric Current I_0 | | $0.805 \pm 0.001 [\mu\text{A}]$ | |
| $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ | $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ |
| 0 | 0.729 | 180 | 0.805 |
| 10 | 0.716 | 190 | 0.779 |
| 20 | 0.668 | 200 | 0.706 |
| 30 | 0.570 | 210 | 0.599 |
| 40 | 0.446 | 220 | 0.470 |
| 50 | 0.316 | 230 | 0.342 |
| 60 | 0.196 | 240 | 0.210 |
| 70 | 0.093 | 250 | 0.101 |
| 80 | 0.029 | 260 | 0.031 |
| 90 | 0.003 | 270 | 0.003 |
| 100 | 0.025 | 280 | 0.023 |
| 110 | 0.091 | 290 | 0.090 |
| 120 | 0.195 | 300 | 0.195 |
| 130 | 0.317 | 310 | 0.310 |
| 140 | 0.460 | 320 | 0.421 |
| 150 | 0.596 | 330 | 0.542 |
| 160 | 0.705 | 340 | 0.649 |
| 170 | 0.779 | 350 | 0.713 |

Table 5: Measurement data for the 1/4-wave plate (rotation angle 0°).

As described in the procedure part, $\sqrt{I/I_0}$ is calculated. Take the first set of data as an example,

$$\sqrt{\frac{I}{I_0}} = \sqrt{\frac{0.729}{0.805}} = 0.9056 \pm 0.0009.$$

Perform similar calculations to each of the other sets of data and the results are presented in Table 6.

| $\theta [^\circ] \pm 2[^\circ]$ | $\sqrt{I/I_0}$ | $\theta [^\circ] \pm 2[^\circ]$ | $\sqrt{I/I_0}$ |
|---------------------------------|---------------------|---------------------------------|---------------------|
| 0 | 0.9516 ± 0.0009 | 180 | 1.0000 ± 0.0009 |
| 10 | 0.9431 ± 0.0009 | 190 | 0.9837 ± 0.0009 |
| 20 | 0.9109 ± 0.0009 | 200 | 0.9365 ± 0.0009 |
| 30 | 0.8415 ± 0.0009 | 210 | 0.8626 ± 0.0009 |
| 40 | 0.7443 ± 0.0010 | 220 | 0.7641 ± 0.0009 |
| 50 | 0.6265 ± 0.0011 | 230 | 0.6518 ± 0.0010 |
| 60 | 0.4934 ± 0.0013 | 240 | 0.5108 ± 0.0013 |
| 70 | 0.3399 ± 0.0018 | 250 | 0.3542 ± 0.0018 |
| 80 | 0.190 ± 0.003 | 260 | 0.196 ± 0.003 |
| 90 | 0.061 ± 0.010 | 270 | 0.061 ± 0.010 |
| 100 | 0.176 ± 0.004 | 280 | 0.169 ± 0.004 |
| 110 | 0.3362 ± 0.0019 | 290 | 0.3344 ± 0.0019 |
| 120 | 0.4922 ± 0.0013 | 300 | 0.4922 ± 0.0013 |
| 130 | 0.6275 ± 0.0011 | 310 | 0.6206 ± 0.0011 |
| 140 | 0.7559 ± 0.0009 | 320 | 0.7232 ± 0.0010 |
| 150 | 0.8604 ± 0.0009 | 330 | 0.8205 ± 0.0009 |
| 160 | 0.9358 ± 0.0009 | 340 | 0.8979 ± 0.0009 |
| 170 | 0.9837 ± 0.0009 | 350 | 0.9411 ± 0.0009 |

Table 6: Results for $\sqrt{I/I_0}$ when rotation angle is 0° .

Then the relationship of $\sqrt{I/I_0}$ and θ are plotted in polar coordinate (Figure 5).




Figure 5: $\sqrt{I/I_0}$ vs. θ relation in polar coordinate when rotation angle is 0° .

4.3.2 Rotation Angle: 20°

The measurement data for 20° rotation angle of 1/4-wave plate are shown in Table 7. Note that when filling the data sheet, 0° and 90° are mistaken and the mistake is corrected in the report.

| Rotation angle of 1/4-wave plate: 20° | | | |
|---------------------------------------|---|---------------------------------|---|
| Maximum Electric Current I_0 | | $0.707 \pm 0.001 [\mu\text{A}]$ | |
| $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ | $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ |
| 0 | 0.589 | 180 | 0.641 |
| 10 | 0.646 | 190 | 0.693 |
| 20 | 0.671 | 200 | 0.707 |
| 30 | 0.662 | 210 | 0.690 |
| 40 | 0.612 | 220 | 0.644 |
| 50 | 0.534 | 230 | 0.567 |
| 60 | 0.453 | 240 | 0.470 |
| 70 | 0.356 | 250 | 0.363 |
| 80 | 0.261 | 260 | 0.263 |
| 90 | 0.179 | 270 | 0.184 |
| 100 | 0.123 | 280 | 0.127 |
| 110 | 0.102 | 290 | 0.104 |
| 120 | 0.116 | 300 | 0.115 |
| 130 | 0.164 | 310 | 0.159 |
| 140 | 0.240 | 320 | 0.224 |
| 150 | 0.345 | 330 | 0.318 |
| 160 | 0.450 | 340 | 0.420 |
| 170 | 0.550 | 350 | 0.513 |

Table 7: Measurement data for the 1/4-wave plate (rotation angle 20°).

Similar to the previous section, $\sqrt{I/I_0}$ is calculated and the results are presented in Table 8 (For sample calculation please refer to section 4.3.1, which will not be repeated here).

| $\theta [^\circ] \pm 2[^\circ]$ | $\sqrt{I/I_0}$ | $\theta [^\circ] \pm 2[^\circ]$ | $\sqrt{I/I_0}$ |
|---------------------------------|---------------------|---------------------------------|---------------------|
| 0 | 0.9127 ± 0.0010 | 180 | 0.9522 ± 0.0010 |
| 10 | 0.9559 ± 0.0010 | 190 | 0.9900 ± 0.0010 |
| 20 | 0.9742 ± 0.0010 | 200 | 1.0000 ± 0.0010 |
| 30 | 0.9677 ± 0.0010 | 210 | 0.9879 ± 0.0010 |
| 40 | 0.9304 ± 0.0010 | 220 | 0.9544 ± 0.0010 |
| 50 | 0.8691 ± 0.0010 | 230 | 0.8955 ± 0.0010 |
| 60 | 0.8005 ± 0.0010 | 240 | 0.8153 ± 0.0010 |
| 70 | 0.7096 ± 0.0011 | 250 | 0.7165 ± 0.0011 |
| 80 | 0.6076 ± 0.0012 | 260 | 0.6099 ± 0.0012 |
| 90 | 0.5032 ± 0.0014 | 270 | 0.5102 ± 0.0014 |
| 100 | 0.4171 ± 0.0017 | 280 | 0.4238 ± 0.0017 |
| 110 | 0.3798 ± 0.0019 | 290 | 0.3835 ± 0.0019 |
| 120 | 0.4051 ± 0.0018 | 300 | 0.4033 ± 0.0018 |
| 130 | 0.4816 ± 0.0015 | 310 | 0.4742 ± 0.0015 |
| 140 | 0.5826 ± 0.0013 | 320 | 0.5629 ± 0.0013 |
| 150 | 0.6986 ± 0.0011 | 330 | 0.6707 ± 0.0012 |
| 160 | 0.7978 ± 0.0011 | 340 | 0.7708 ± 0.0011 |
| 170 | 0.8820 ± 0.0010 | 350 | 0.8518 ± 0.0010 |

Table 8: Results for $\sqrt{I/I_0}$ when rotation angle is 20° .

Then the relationship of $\sqrt{I/I_0}$ and θ are plotted in polar coordinate (Figure 6).




Figure 6: $\sqrt{I/I_0}$ vs. θ relation in polar coordinate when rotation angle is 20° , together with the maximum point of rotation angle 70° (marked in red).

4.3.3 Rotation Angle: 45°

The measurement data for 45° rotation angle of 1/4-wave plate are shown in Table 9. Note that when filling the data sheet, 0° and 90° are mistaken and the mistake is corrected in the report.

| Rotation angle of 1/4-wave plate: 45° | | | |
|---------------------------------------|---|---------------------------------|---|
| Maximum Electric Current I_0 | | $0.395 \pm 0.001 [\mu\text{A}]$ | |
| $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ | $\theta [^\circ] \pm 2[^\circ]$ | $I [\mu\text{A}] \pm 0.001 [\mu\text{A}]$ |
| 0 | 0.359 | 180 | 0.392 |
| 10 | 0.360 | 190 | 0.387 |
| 20 | 0.356 | 200 | 0.378 |
| 30 | 0.354 | 210 | 0.370 |
| 40 | 0.351 | 220 | 0.368 |
| 50 | 0.348 | 230 | 0.365 |
| 60 | 0.354 | 240 | 0.363 |
| 70 | 0.361 | 250 | 0.365 |
| 80 | 0.364 | 260 | 0.369 |
| 90 | 0.371 | 270 | 0.375 |
| 100 | 0.375 | 280 | 0.386 |
| 110 | 0.381 | 290 | 0.392 |
| 120 | 0.385 | 300 | 0.389 |
| 130 | 0.390 | 310 | 0.378 |
| 140 | 0.393 | 320 | 0.366 |
| 150 | 0.395 | 330 | 0.365 |
| 160 | 0.393 | 340 | 0.369 |
| 170 | 0.395 | 350 | 0.364 |

Table 9: Measurement data for the 1/4-wave plate (rotation angle 45°).

Similar to the previous section, $\sqrt{I/I_0}$ is calculated and the results are presented in Table 10 (For sample calculation please refer to section 4.3.1, which will not be repeated here).

| $\theta [^\circ] \pm 2[^\circ]$ | $\sqrt{I/I_0}$ | $\theta [^\circ] \pm 2[^\circ]$ | $\sqrt{I/I_0}$ |
|---------------------------------|---------------------|---------------------------------|---------------------|
| 0 | 0.9533 ± 0.0018 | 180 | 0.9962 ± 0.0018 |
| 10 | 0.9547 ± 0.0018 | 190 | 0.9898 ± 0.0018 |
| 20 | 0.9494 ± 0.0018 | 200 | 0.9782 ± 0.0018 |
| 30 | 0.9467 ± 0.0018 | 210 | 0.9678 ± 0.0018 |
| 40 | 0.9427 ± 0.0018 | 220 | 0.9652 ± 0.0018 |
| 50 | 0.9386 ± 0.0018 | 230 | 0.9613 ± 0.0018 |
| 60 | 0.9467 ± 0.0018 | 240 | 0.9586 ± 0.0018 |
| 70 | 0.9560 ± 0.0018 | 250 | 0.9613 ± 0.0018 |
| 80 | 0.9600 ± 0.0018 | 260 | 0.9665 ± 0.0018 |
| 90 | 0.9691 ± 0.0018 | 270 | 0.9744 ± 0.0018 |
| 100 | 0.9744 ± 0.0018 | 280 | 0.9885 ± 0.0018 |
| 110 | 0.9821 ± 0.0018 | 290 | 0.9962 ± 0.0018 |
| 120 | 0.9873 ± 0.0018 | 300 | 0.9924 ± 0.0018 |
| 130 | 0.9937 ± 0.0018 | 310 | 0.9782 ± 0.0018 |
| 140 | 0.9975 ± 0.0018 | 320 | 0.9626 ± 0.0018 |
| 150 | 1.0000 ± 0.0018 | 330 | 0.9613 ± 0.0018 |
| 160 | 0.9975 ± 0.0018 | 340 | 0.9665 ± 0.0018 |
| 170 | 1.0000 ± 0.0018 | 350 | 0.9600 ± 0.0018 |

Table 10: Results for $\sqrt{I/I_0}$ when rotation angle is 45° .

Then the relationship of $\sqrt{I/I_0}$ and θ are plotted in polar coordinate (Figure 7).




Figure 7: $\sqrt{I/I_0}$ vs. θ relation in polar coordinate when rotation angle is 45° .

To compare the result with circular polarization, as described in the procedure part, the data is also plotted in Cartesian coordinate and linear fit is performed (Figure 8). The slope of the linear fitting is $6 \times 10^{-5} \pm 6 \times 10^{-5}$.




Figure 8: $\sqrt{I/I_0}$ vs. θ relation in the Cartesian coordinate.

4.3.4 Rotation Angle: 70°

In the condition of the rotation angle is 70°, the data when the light intensity reaches its maximum are presented in Table 11. Still, the data is subtracted by 90°. The position of this point is marked in red in Figure 6. It can be seen that the value of θ when the light intensity reaches maximum in the case of 70° rotation angle is nearly symmetric about $\theta = 0$ axis with that ($\theta = 20^\circ(200^\circ)$) in the case of 20° rotation angle.

| Rotation angle of the 1/4-wave plate: 70° | |
|---|-------|
| $\theta [^\circ] \pm 2[^\circ]$ | 161 |
| $I [\mu A] \pm 0.001 [\mu A]$ | 0.682 |

Table 11: Measurement data for the 1/4-wave plate (rotation angle 70°).

5 Discussion and Conclusions

5.1 Demonstration of Malus' Law

The slope of the fitting is 1.01 ± 0.02 and the Pearson's r is 0.99916, which is very close to 1. This suggests that the value of I/I_0 is proportional to the value of $\cos^2 \theta$ with the coefficient to be about 1.

Theoretically, by Eq. (1), the slope is

$$\frac{I/I_0}{\cos^2 \theta} = 1,$$

Therefore, the relative error is

$$\varepsilon = \frac{1.01 - 1}{1} \times 100\% = 1\%.$$

This in an acceptable range of error verifies the Malus' Law.

5.2 Linearly Polarized Light and the Half-wave Plate

The slope of the linear fit is 2.01 ± 0.02 . Theoretically, as introduced in the Introduction part, the rotation angle of the polarization axis is twice of the origin angle for a half-wave plate. Therefore the theoretical value of the slope of our linear fitting is 2. The relative error is therefore

$$\varepsilon = \frac{2.01 - 2}{2} \times 100\% = 0.5\%.$$

This conforms to the theoretical fact.

5.3 Circularly and Elliptically Polarized Light and the 1/4-wave Plate

5.3.1 Rotation Angle: 0°

From Table 5 and Figure 5, it can be seen that the maximum of light intensity occurs at about $\theta = 0^\circ$. This suggests that polarizing axis is parallel to the optical axis of the plate, which is the theoretical conclusion stated in the Introduction part. The shape of the plot also indicates that it is linearly polarized.

5.3.2 Rotation Angle: 20°

From Table 7 and Figure 6, it can be seen that the maximum of light intensity occurs at about $\theta = 20^\circ$. This suggests that polarizing axis forms 20° to the optical axis of the plate, which conforms to the fact stated in [2]. The shape of the plot also indicates that it is elliptically polarized.

5.3.3 Rotation Angle: 45°

From Figure 7, it can be seen that the value of $\sqrt{I/I_0}$ changes little when θ changes. Besides, the slope of linear fit of the plot of $\sqrt{I/I_0}$ vs. θ in the Cartesian coordinate in Figure 8 is $6 \times 10^{-5} \pm 6 \times 10^{-5}$ and the Pearson's coefficient is only 0.328. All the results suggest that $\sqrt{I/I_0}$ is a constant, with no relation with the value of θ . Therefore it can be concluded that when rotation angle is 45° , the light is circularly polarized. This conforms to the fact stated in the Introduction part.

5.3.4 Rotation Angle: 70°

From section 4.3.4, we can state that the value of θ when the light intensity reaches maximum in the case of 70° rotation angle is nearly symmetric about $\theta = 0$ axis with that ($\theta = 20^\circ(200^\circ)$) in the case of 20° rotation angle. This conclusion conforms to the theoretical conclusion stated in [2].

5.4 Causes for errors and uncertainties

Possible causes for errors and uncertainties are listed below:

1. The light of our torch may affect the light intensity detected by the photo-cell.
2. The precision of the scale on the optical elements is 2° , which results in a relatively large uncertainty.

3. It is really tough to adjust the center of all the optical elements to be in a line and we failed to do that in our experiment. This may lead to errors to the experiment.

5.5 Suggestions

Some suggestions for further improvement:

1. The experiment may be performed in an environment that is absolutely dark, where the readings are also available.
2. Improve the precision of the scale on the optical elements.

5.6 Conclusions

In this lab, we studied the polarization phenomenon and verify Malus' law. The way half- and quarter-wave plates work in optical systems is also explored. The experimental results we obtained is of relatively small uncertainties and conform to the theoretical facts in an acceptable range of uncertainty.

References

- [1] VP241 Exercise 4: Polarization of Light, Department of Physics, Shanghai Jiaotong University.
- [2] University Physics, Section 33.5, Young and Freedman.

A Measurement Uncertainty Analysis

A.1 Uncertainty for Data in Demonstration of Malus' Law

The uncertainty of $\cos^2 \theta$ is

$$\begin{aligned} u_{\cos^2 \theta} &= \sqrt{\left(\frac{\partial \cos^2 \theta}{\partial \theta} u_\theta\right)^2} \\ &= |-2 \cos \theta \sin \theta u_\theta| \\ &= |\sin 2\theta u_\theta|, \end{aligned}$$

where $u_\theta = 2[^\circ] = \pi/90$. Take $\theta = 0$ as an example,

$$u_{\cos^2 \theta} = |\sin 2\theta u_\theta| = |\sin(2 \times 5^\circ) \times \frac{\pi}{90}| = 0.006,$$

All the results are shown in Table 12.

The uncertainty of I/I_0 is

$$\begin{aligned} u_{I/I_0} &= \sqrt{\left(\frac{\partial I/I_0}{\partial I} u_I\right)^2 + \left(\frac{\partial I/I_0}{\partial I_0} u_{I_0}\right)^2} \\ &= \sqrt{\left(\frac{u_I}{I_0}\right)^2 + \left(-\frac{I}{I_0^2} u_{I_0}\right)^2}, \end{aligned}$$

where $u_I = u_{I_0} = 0.001 [\mu\text{A}]$, $I_0 = 1.037 \pm 0.001 [\mu\text{A}]$. Take the first set of data as an example,

$$u_{I/I_0} = \sqrt{\left(\frac{u_I}{I_0}\right)^2 + \left(-\frac{I}{I_0^2} u_{I_0}\right)^2} = \sqrt{\left(\frac{0.001}{1.037}\right)^2 + \left(-\frac{1.037}{1.037^2} \times 0.001\right)^2} = 0.0013,$$

All the results are shown in Table 12.

| $u_{\cos^2 \theta}$ | u_{I/I_0} | $u_{\cos^2 \theta}$ | u_{I/I_0} |
|---------------------|-------------|---------------------|-------------|
| 0 | 0.0013 | 0.03 | 0.0011 |
| 0.006 | 0.0013 | 0.03 | 0.0010 |
| 0.012 | 0.0013 | 0.03 | 0.0010 |
| 0.017 | 0.0013 | 0.03 | 0.0010 |
| 0.02 | 0.0013 | 0.02 | 0.0010 |
| 0.03 | 0.0013 | 0.017 | 0.0010 |
| 0.03 | 0.0012 | 0.012 | 0.0010 |
| 0.03 | 0.0012 | 0.006 | 0.0010 |
| 0.03 | 0.0011 | 0 | 0.0010 |
| 0.03 | 0.0011 | | |

Table 12: Uncertainty of $\cos^2 \theta$ and I/I_0 .

A.2 Uncertainty for Linearly Polarized Light and the Half-wave Plate

In this part, the uncertainty of the measurement depends on the device, which is

$$u_\theta = 2^\circ.$$

A.3 Uncertainty for Circularly and Elliptically Polarized Light and the 1/4-wave Plate

The uncertainty of $\sqrt{I/I_0}$ is

$$\begin{aligned} u_{\sqrt{I/I_0}} &= \sqrt{\left(\frac{\partial\sqrt{I/I_0}}{\partial I}u_I\right)^2 + \left(\frac{\partial\sqrt{I/I_0}}{\partial I_0}u_{I_0}\right)^2} \\ &= \sqrt{\frac{1}{4II_0}u_I^2 + \frac{I}{4I_0^3}u_{I_0}^2}, \end{aligned}$$

where $u_I = u_{I_0} = 0.001$ [μA], $I_0 = 0.805 \pm 0.001$ [μA] for rotation angle of 0° , $I_0 = 0.707 \pm 0.001$ [μA] for rotation angle of 20° , $I_0 = 0.395 \pm 0.001$ [μA] for rotation angle of 45° . Take the first set of data as an example,

$$u_{\sqrt{I/I_0}} = \sqrt{\frac{1}{4 \times 0.729 \times 0.805} \times 0.001^2 + \frac{0.729}{4 \times 0.805^3} \times 0.001^2} = 0.0009.$$

All the results are shown in Table 13.

| Rotation angle: 0° | Rotation angle: 20° | Rotation angle: 45° |
|--------------------|---------------------|---------------------|
| $u \sqrt{I/I_0}$ | $u \sqrt{I/I_0}$ | $u \sqrt{I/I_0}$ |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0010 | 0.0010 | 0.0018 |
| 0.0011 | 0.0010 | 0.0018 |
| 0.0013 | 0.0010 | 0.0018 |
| 0.0018 | 0.0011 | 0.0018 |
| 0.003 | 0.0012 | 0.0018 |
| 0.010 | 0.0014 | 0.0018 |
| 0.004 | 0.0017 | 0.0018 |
| 0.0019 | 0.0019 | 0.0018 |
| 0.0013 | 0.0018 | 0.0018 |
| 0.0011 | 0.0015 | 0.0018 |
| 0.0009 | 0.0013 | 0.0018 |
| 0.0009 | 0.0011 | 0.0018 |
| 0.0009 | 0.0011 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |
| 0.0010 | 0.0010 | 0.0018 |
| 0.0013 | 0.0010 | 0.0018 |
| 0.0018 | 0.0011 | 0.0018 |
| 0.003 | 0.0012 | 0.0018 |
| 0.010 | 0.0014 | 0.0018 |
| 0.004 | 0.0017 | 0.0018 |
| 0.0019 | 0.0019 | 0.0018 |
| 0.0013 | 0.0018 | 0.0018 |
| 0.0011 | 0.0015 | 0.0018 |
| 0.0010 | 0.0013 | 0.0018 |
| 0.0009 | 0.0012 | 0.0018 |
| 0.0009 | 0.0011 | 0.0018 |
| 0.0009 | 0.0010 | 0.0018 |

Table 13: Uncertainty for $\sqrt{I/I_0}$ when the rotation angle is 0°, 20° and 45°.

B Data Sheet

Please find the original data sheet attached at the end of the report.

UM-SJTU PHYSICS LABORATORY
DATA SHEET (EXERCISE 4)

Name: 康家铭

Student ID: 518021911220

Name: 洪锡辰 (Hong Xichen)

Student ID: 518373910011

Group: 17

Date: 8. Nov. 2019

NOTICE. Please remember to show the data sheet to your instructor before leaving the laboratory. The data sheet will not be accepted if the data are recorded with pencil or modified by correction fluid/tape. If a mistake is made in recording a datum item, cancel the wrong value by drawing a fine line through it, record the correct value legibly, and ask your instructor to confirm the correction. Please remember to take a record of the precision of the instruments used. You are required to hand in the original data with your lab report, so please keep the data sheet properly.

1.418

Uncertainty of θ is $[2]^\circ$.

| Maximum Electric Current I_0 | | | |
|--------------------------------|---------------------------|-------------------------|---------------------------|
| θ | $1 \mu A \pm 0.001 \mu A$ | $1.037 \pm 0.001 \mu A$ | $1 \mu A \pm 0.001 \mu A$ |
| 0° | 1.037 | 50° | 0.470 |
| 5° | 1.033 | 55° | 0.373 |
| 10° | 1.017 | 60° | 0.294 |
| 15° | 0.986 | 65° | 0.226 |
| 20° | 0.947 | 70° | 0.147 |
| 25° | 0.896 | 75° | 0.085 |
| 30° | 0.852 | 80° | 0.042 |
| 35° | 0.730 | 85° | 0.015 |
| 40° | 0.643 | 90° | 0.003 |
| 45° | 0.560 | | |

Table 1. Measurement data Malus' law demonstration.

Instructor's signature: EXX

| Rotation angle of the 1/2-wave plate | Rotation angle of the analyzer [$^{\circ}$] \pm <u>2</u> $^{\circ}$ |
|--------------------------------------|---|
| initial | 210 |
| 10 $^{\circ}$ | 232 |
| 20 $^{\circ}$ | 250 |
| 30 $^{\circ}$ | 271 |
| 40 $^{\circ}$ | 292 |
| 50 $^{\circ}$ | 312 |
| 60 $^{\circ}$ | 331 |
| 70 $^{\circ}$ | 352 |
| 80 $^{\circ}$ | 12 |
| 90 $^{\circ}$ | 36 |

Table 2. Measurement data for the 1/2-wave plate.

Instructor's signature: 8xx

maximum

| Rotation angle of 1/4-wave plate: 0° | | | |
|--------------------------------------|---------------------------|-------------------------|---------------------------|
| Maximum Electric Current I_0 | | $0.805 \pm 0.001 \mu A$ | |
| θ | $I \mu A \pm 0.001 \mu A$ | θ | $I \mu A \pm 0.001 \mu A$ |
| 0° | 0.003 | 180° | 0.003 |
| 10° | 0.023 | 190° | 0.025 |
| 20° | 0.090 | 200° | 0.091 |
| 30° | 0.195 | 210° | 0.195 |
| 40° | 0.310 | 220° | 0.317 |
| 50° | 0.421 | 230° | 0.460 |
| 60° | 0.542 | 240° | 0.596 |
| 70° | 0.649 | 250° | 0.705 |
| 80° | 0.713 | 260° | 0.779 |
| 90° | 0.729 | 270° | 0.805 |
| 100° | 0.716 | 280° | 0.779 |
| 110° | 0.668 | 290° | 0.706 |
| 120° | 0.570 | 300° | 0.599 |
| 130° | 0.446 | 310° | 0.470 |
| 140° | 0.316 | 320° | 0.342 |
| 150° | 0.196 | 330° | 0.210 |
| 160° | 0.093 | 340° | 0.101 |
| 170° | 0.029 | 350° | 0.031 |

Table 3. Measurement data for the 1/4-wave plate (rotation angle 0°).

Table 3-5, 0° to 180°

Instructor's signature: 87X

Sn 454 20°

| Rotation angle of the 1/4-wave plate: 20° | | | |
|---|---------------------------|-------------------------|---------------------------|
| Maximum Electric Current I_0 | | $0.101 \pm 0.001 \mu A$ | |
| θ | $I \mu A \pm 0.001 \mu A$ | θ | $I \mu A \pm 0.001 \mu A$ |
| 0° | 0.184 | 180° | 0.179 |
| 10° | 0.127 | 190° | 0.123 |
| 20° | 0.100 | 200° | 0.102 |
| 30° | 0.115 | 210° | 0.116 |
| 40° | 0.159 | 220° | 0.164 |
| 50° | 0.224 | 230° | 0.240 |
| 60° | 0.318 | 240° | 0.345 |
| 70° | 0.418 0.420 | 250° | 0.450 |
| 80° | 0.513 | 260° | 0.550 |
| 90° | 0.589 | 270° | 0.641 |
| 100° | 0.646 | 280° | 0.673 |
| 110° | 0.671 | 290° | 0.707 |
| 120° | 0.662 | 300° | 0.690 |
| 130° | 0.612 | 310° | 0.644 |
| 140° | 0.534 | 320° | 0.567 |
| 150° | 0.453 | 330° | 0.470 |
| 160° | 0.356 | 340° | 0.363 |
| 170° | 0.261 | 350° | 0.265 0.263 |

Table 4. Measurement data for the 1/4-wave plate (rotation angle 20°).

Instructor's signature: XXX

| Rotation angle of the 1/4-wave plate: 45° | | | |
|---|-------------------------------|-------------------------|-------------------------------|
| Maximum Electric Current I_0 | | $0.375 \pm 0.001 \mu A$ | |
| θ | $I [\mu A] \pm 0.001 [\mu A]$ | θ | $I [\mu A] \pm 0.001 [\mu A]$ |
| 0° | 0.375 | 180° | 0.371 |
| 10° | 0.386 | 190° | 0.375 |
| 20° | 0.392 | 200° | 0.381 |
| 30° | 0.389 | 210° | 0.385 |
| 40° | 0.378 | 220° | 0.390 |
| 50° | 0.366 | 230° | 0.393 |
| 60° | 0.365 | 240° | 0.395 |
| 70° | 0.369 | 250° | 0.393 |
| 80° | 0.359 | 260° | 0.395 |
| 90° | 0.360 | 270° | 0.392 |
| 100° | 0.356 | 280° | 0.387 |
| 110° | 0.354 | 290° | 0.378 |
| 120° | 0.354 | 300° | 0.370 |
| 130° | 0.351 | 310° | 0.368 |
| 140° | 0.348 | 320° | 0.365 |
| 150° | 0.354 | 330° | 0.363 |
| 160° | 0.361 | 340° | 0.365 |
| 170° | 0.360 | 350° | 0.369 |

Table 5. Measurement data for the 1/4-wave plate (rotation angle 45°).

| Rotation angle of the 1/4-wave plate: 70° | |
|---|-------|
| $\theta [^\circ] \pm [2]^\circ$ | 251° |
| $I [\mu A] \pm 0.001 [\mu A]$ | 0.682 |

initial 29°
280°

Table 6. Measurement data for the 1/4-wave plate (rotation angle 70°).

Instructor's signature: _____

8xx